

# **AE-1250 – Roll Up Intubation Cases Cleaning Recommendations**

The cases that are made from the cloth material should withstand machine washing with no problems. Several customers have cleaned their products in a washing machine and have not reported any difficulties. I believe they generally let them air dry rather than using a clothes dryer. We can't make any representation as to the infection control / body fluid problem, because we have no way of monitoring what procedures the end user employs.

Be aware that on products that are made from the antibacterial / vinyl material, the user can spray the vinyl with Lysol, 409, etc., but the bags are sewn together with cloth thread. This thread will absorb body fluids and any other liquid. Obviously, the elastic material on the inside will also absorb liquids. The bottom line is that there is no such thing as a bag or case that will repel liquids 100%.

In summary, it's a good idea to spot clean or wash the product, but there is no guarantee that the cleaning will make it germ free. From a sales point of view we encourage the customers to burn straps and bags if they have concerns about contagion.